



aCalendar	Simple calendar
theWiz	This is a mini browser that understands special MaxTRAQ hyperlinks that can be used to control MaxTRAQ from theWiz. The intention of theWiz is to provide a simple way to create wizards or student course material that the user can follow along. The default file extension for theWiz is HTM or HTML. See the scripting tutorial for more info.
aCalculator	Scientific calculator
myNotes	A notepad to allow you to keep notes. The notes are saved as text files. You can drag and drop a text file onto MaxTRAQ to open up the file in myNotes.
eventREC	Manual Event Recorder. Allows you to step through the video and record events such as heel down/toe off etc.

- ◆ Teachers can use their own course material in theWiz
- ◆ MaxTRAQ a great Mocap tool to complement current curriculum
- ◆ MaxTRAQ is used worldwide as a teaching tool/lab course